


Results and conclusion²

According to **DIN EN 17126 (2019)** and the “Requirements and methods for **VAH**-certification of chemical disinfection processes” (**2015/2025**), the batches B25002/ B25001 of the 2-component-product **MEIKO® ACTIVE SPOROSAN® 8100 A+B**, when applied at the concentration / contact time - relations of at least **80 % / 15 s** under clean conditions (0.03 % albumin), **possesses sporicidal efficacy** ($\log_{10} \text{RF} \geq 4$) for the reference strains *B. subtilis* and *B. cereus* (Tab. 1.1 – 2.2).

According to **DIN EN 17126 (2019)** and the “Requirements and methods for **VAH**-certification of chemical disinfection processes” (**2015/2025**), the batches B25002/ B25001 of the 2-component-product **MEIKO® ACTIVE SPOROSAN® 8100 A+B**, when applied at the concentration / contact time - relations of at least **80 % / 15 s** under clean conditions (0.03 % albumin), **possesses sporicidal efficacy against *C. difficile*** ($\log_{10} \text{RF} \geq 4$) for the reference strain *C. difficile* (Tab. 3.1 – 3.2).

Results are considered validated in accordance with the requirements of the **DIN EN 17126 (2019)** and **VAH (2015/2025)**.

Greifswald, April 01, 2026


Dr. rer. med. (Dipl. Biol.) T. Koburger-Janssen
- General Manager -


Prof. Dr. med. A. Kramer
- MD for Hygiene and Environmental Medicine -