


Results and conclusion ²

In accordance with the **EN 17387 (2021)**, the batches H 023005/2 and N 023005/1 of the 2-component-disinfectant **SPOROSAN® WD 42**, when applied as described above at the concentration / contact time - relation of at least **100 % / 1 min** at $20 \pm 1 ^\circ\text{C}$ under clean conditions (0.03 % albumin), **possesses sporicidal efficacy** ($\log_{10}$ RF ≥ 4) for reference strains *B. subtilis*, *B. cereus* and *C. difficile* (Tab. 1.1 - 3).

Results are considered validated in accordance with the requirements of **EN 17387 (2021)**.

Greifswald, April 26, 2024


Dr. rer. med. (Dipl. Biol.) T. Koburger-Janssen
- General Manager -


Prof. Dr. med. A. Kramer
MD for Hygiene and Environmental Medicine -

